

Resolvigen4: summary of changes from 0.0.1.7 to 0.0.2.0

Written by: Carlo Magnaghi

Date : 25 May 2017

Errors fixed

0.0.1.7_S1

Antigenic expressions leading to the exclusion of an anti-k (anti-Cellano) antibody are now displayed correctly.

0.0.1.7_S2

The Vision/Innova import module no longer loses the link to the Panels expiration date when editing links. This caused Resolvigen to ask panel information again to the user.

0.0.1.7_S3

It is now possible to edit donor's RBCs reactions when editing data inside a test.

0.0.1.7_S4

The option to prevent exporting personal data when exporting test data now works correctly.

0.0.1.7_S5

When removing a patient from the archives, sometimes Resolvigen offered the user the option to delete all tests instead of only the tests for the given patient. This is now fixed.

0.0.1.7_S6

Sometimes, when in presence of incoherent scores, so that Resolvigen decided to proceed with the interpretation ignoring score differences, Resolvigen gave anyway a suggestion of antibodies best matching strong scores, which ended up to be misleading. This has been fixed.

0.0.1.7_S7

Now Resolvigen shows exclusion reasons even for scores given only in the IgG phase.

0.0.1.7_S8

When merging results from Innova sample-id was not displayed; this is now fixed.

0.0.1.7_S9

The list of rare donor's RBCs is now always updated immediately after changing filter criteria.

Change requests implemented

0.0.1.7_R2

When giving a result obtained with Fisher formula Resolvigen displayed the correlated phenotype expression with + and - (e.g. Jka-); the '-' is hardly visible; now Resolvigen displays 'positive' and 'negative' instead.

0.0.1.7_R3

Some minor fixes in the German translation of the software

0.0.1.7_R4

When displaying antibodies in the procedures for the identification against public antigens, now hyperlinks to the online documentation are shown also in the preliminary page: its used to be only in the final page.

0.0.1.7_R5

When displaying the antigram for patient's and donor's RBCs now the ABO phenotype is also shown.

0.0.1.7_R6

Now it is possible to use custom Panels with a number of RBCs from 2 to 15.

A code must be supplied to customers to enable them to edit custom panels.

Given the large number of panel types now possible, now when selecting a panel type for adding a panel to a test, only the panels types actually present in the archives are displayed.

0.0.1.7_R7

When asking the user to merge results, now YES/NO buttons are shown instead of OK/Cancel.

NOTE: some improvements have also been implemented in the diagnostic engine.

WARNING: Resolvigen 4 0.0.2.0 is compatible with Vision software versions up to 4.8; validation with versions 5.x is pending.